

Streampix and Troublepix: Digital Video Recording Software and Solution



*Fast and easy to use
Direct to your laptop*

*Portability
and long
recording periods*

*Compatible
with many cameras*

*Up to 200 frames
per seconds
at 640 x 480*

*Laptop with GigE
or Firewire B camera
running at 200 fps.*

*Check with NorPix
for suitable laptop
configuration.*

*Turnkey solutions
available.*

NorPix, Inc.

Recording to RAM

<i>Resolution</i>	<i>Frames per second</i>	<i>4 Gbytes</i>	<i>3 Gbytes</i>	<i>2 Gbytes</i>
640 x 480	200	45 s	28 s	12 s
320 x 240	371	1 min 44 s	1 min 6 s	28 s
240 x 160	512	2 min 22 s	1 min 30 s	28 s

Direct to disk

<i>Resolution</i>	<i>Frames per second</i>	<i>Medium</i>	<i>Length of Replay</i>
640 x 480	200	to hard disk	20 min
320 x 240	371	to hard disk	1 h
240 x 160	512	to hard disk	1 h 47 min

Using JPEG compression with 75% quality

<i>Resolution</i>	<i>Frames per second</i>	<i>Medium</i>	<i>Recording to RAM</i>
640 x 480	200	16 h to hard disk	20 min

NorPix, Inc.
1751 Richardson, suite 6117
Montreal, Canada
H3K 1G6
sales@norpix.com

Tel: (514) 846- 0009 x 23
Fax: (514) 846- 0117
<http://www.norpix.com>